

Antibodies Product Datasheet

Product Name:	FGF-1 Antibody
Catalog Number:	R12-2141
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FGF-1. AA range:7-56
Concentration:	1 mg/ml
Gene Name:	FGF1
Protein Name:	Fibroblast growth factor 1
Other Name:	FGF1; FGFA; Fibroblast growth factor 1; FGF-1; Acidic fibroblast growth factor; aFGF; Endothelial cell growth factor; ECGF; Heparin-binding growth factor 1; HBGF-1
Human Swiss Prot No.:	P05230
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FGF-1 Polyclonal Antibody detects endogenous levels of FGF-1 protein.
Storage and Stability:	-20°C/1 year