

## Antibodies Product Datasheet

<b>Product Name:</b>	Synuclein- $\alpha$ Antibody
<b>Catalog Number:</b>	C30743
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from the Internal region of human SNCA. AA range:21-70
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SNCA
<b>Protein Name:</b>	Alpha-synuclein
<b>Other Name:</b>	SNCA; NACP; PARK1; Alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; NACP
<b>Human Swiss Prot No.:</b>	P37840
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Synuclein- $\alpha$ Polyclonal Antibody detects endogenous levels of Synuclein- $\alpha$ protein.
<b>Storage and Stability:</b>	-20°C/1 year