

Antibodies Product Datasheet

Product Name:	p56 Dok-2 Antibody
Catalog Number:	R12-3290
Reactivity:	Human, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	p56 Dok-2 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human p56 Dok-2.
Concentration:	1 mg/ml
Gene Name:	DOK2
Other Name:	Docking protein 2; Dok-R; Dok-related protein; DOK2; DokR; downstream of tyrosine kinase 2; FRIP; IL-4 receptor interacting protein; p56(dok-2)
Human Swiss Prot No.:	O60496
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	p56 Dok-2 Antibody antibody detects endogenous levels of p56 Dok-2.
Storage and Stability:	Stable at -20°C for at least 1 year.