

Antibodies Product Datasheet

Product Name:	AP3B2 Antibody
Catalog Number:	C49450
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AP3B2
Concentration:	1 mg/ml
Gene Name:	AP3B2
Protein Name:	AP-3 complex subunit beta-2 (Adapter-related protein complex 3 subunit beta-2) (Adaptor protein complex AP-3 subunit beta-2) (Beta-3B-adaptin) (Clathrin assembly protein complex 3 beta-2 large chain) (Neuron-specific vesicle coat protein beta-NAP)
Other Name:	
Human Swiss Prot No.:	Q13367
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of AP3B2 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year