
Antibodies Product Datasheet

Product Name:	c-Myb (Acetyl-Lys480) Antibody
Catalog Number:	D12255
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human c-Myb (Acetyl Lys480)
Concentration:	1 mg/ml
Gene Name:	MYB
Protein Name:	c-Myb (Acetyl Lys480)
Other Name:	Transcriptional activator Myb (Proto-oncogene c-Myb)
Human Swiss Prot No.:	P10242
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse c-Myb (Acetyl Lys480)
Storage and Stability:	-20°C/1 year