

Antibodies Product Datasheet

Product Name:	FADS2 Antibody
Catalog Number:	C49043
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FADS2
Concentration:	1 mg/ml
Gene Name:	FADS2
Protein Name:	Fatty acid desaturase 2 (EC 1.14.19.-) (Delta(6) fatty acid desaturase) (D6D) (Delta(6) desaturase) (Delta-6 desaturase)
Other Name:	
Human Swiss Prot No.:	O95864
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FADS2 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year