

Antibodies Product Datasheet

Product Name:	SETD1A Antibody
Catalog Number:	C51037
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SETD1A
Concentration:	1 mg/ml
Gene Name:	SETD1A KIAA0339 KMT2F SET1 SET1A
Protein Name:	Histone-lysine N-methyltransferase SETD1A (EC 2.1.1.43) (Lysine N-methyltransferase 2F) (SET domain-containing protein 1A) (hSET1A) (Set1/Ash2 histone methyltransferase complex subunit SET1)
Other Name:	
Human Swiss Prot No.:	O15047
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SETD1A at Human
Storage and Stability:	Store at -20° C/1 year