

## Antibodies Product Datasheet

<b>Product Name:</b>	DEN5A Antibody
<b>Catalog Number:</b>	C43421
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human protein . at AA range: 650-730
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	DENND5A KIAA1091 RAB6IP1
<b>Protein Name:</b>	DENN domain-containing protein 5A (Rab6-interacting protein 1) (Rab6IP1)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q6IQ26
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Specificity:</b>	DEN5A Polyclonal Antibody detects endogenous levels of protein.
<b>Storage and Stability:</b>	-20°C/1 year