

Antibodies Product Datasheet

Product Name:	CDK5 Mouse mAb(3D1)
Catalog Number:	N1439
Reactivity:	Human, Mouse,Rat
Applications:	IHC-p,WB
Dilutions:	IHC-p1:50-200 ,WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human CDK5
Concentration:	1 mg/ml
Gene Name:	CDK5 CDKN5
Protein Name:	CDK5
Other Name:	Cyclin-dependent kinase 5 (EC 2.7.11.22) (Cell division protein kinase 5) (Serine/threonine-protein kinase PSSALRE) (Tau protein kinase II catalytic subunit) (TPKII catalytic subunit)
Human Swiss Prot No.:	Q00535
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.04% sodium azide.
Specificity:	This antibody detects endogenous levels of CDK5 at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year