

Antibodies Product Datasheet

Product Name:	NPRL2 Antibody
Catalog Number:	C20645
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human NPRL2.
Concentration:	1 mg/ml
Gene Name:	NPRL2
Protein Name:	Nitrogen permease regulator 2-like protein
Other Name:	NPRL2; TUSC4; Nitrogen permease regulator 2-like protein; NPR2-like protein; Gene 21 protein; G21 protein; Tumor suppressor candidate 4
Human Swiss Prot No.:	Q8WTW4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NPRL2 Polyclonal Antibody detects endogenous levels of NPRL2 protein.
Storage and Stability:	-20°C/1 year