

Antibodies Product Datasheet

Product Name:	Anti-Human Cytokeratin 13 (ABT151) Mouse Monoclonal Antibody
Catalog Number:	V0286
Clone ID:	ABT151
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	KRT13
Gene Name:	P13646
Other Name:	47 kDa cytokeratin, CK-13, CK13, Cytokeratin 13, Cytokeratin-13, K13, K1C13_HUMAN, Ka13, Keratin 13, Keratin, keratin type I cytoskeletal 13, Keratin-13, Krt-1.13, Krt1-13, KRT13, MGC161462, MGC3781, type I cytoskeletal 13, Type I keratin Ka13, WSN2
Isotype	IgG1, kappa
Localization:	Tonsil
Positive Tissue Control:	Cytoplasmic, Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human CK13 protein. In western blotting of A431 cell lysate, the antibody can label a 50 kDa band corresponding to CK13.
Storage and Stability:	-20°C/1 year