

Antibodies Product Datasheet

Product Name:	ARX Antibody
Catalog Number:	C20664
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ARX . at AA range: 250-330
Concentration:	1 mg/ml
Gene Name:	ARX
Protein Name:	Homeobox protein ARX
Other Name:	ARX; Homeobox protein ARX; Aristaless-related homeobox
Human Swiss Prot No.:	Q96QS3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARX Polyclonal Antibody detects endogenous levels of ARX protein.
Storage and Stability:	-20°C/1 year