

Antibodies Product Datasheet

Product Name:	AhR Antibody
Catalog Number:	R12-2466
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	AhR Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AhR.
Concentration:	1 mg/ml
Gene Name:	AHR
Other Name:	Aryl hydrocarbon receptor; Ah receptor; AhR; Class E basic helix-loop-helix protein 76; bHLHe76; BHLHE76
Human Swiss Prot No.:	P35869/A9YTQ3
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	AhR Antibodyantibody detects endogenous levels of AhR.
Storage and Stability:	Stable at -20°C for at least 1 year.