

Antibodies Product Datasheet

Product Name:	TGase2 Antibody
Catalog Number:	C0350
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Transglutaminase 2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	TGM2
Protein Name:	Protein-glutamine gamma-glutamyltransferase 2
Other Name:	TGM2; Protein-glutamine gamma-glutamyltransferase 2; Tissue transglutaminase; Transglutaminase C; TG(C); TGC; TGase C; Transglutaminase H; TGase H; Transglutaminase-2; TGase-2
Human Swiss Prot No.:	P21980
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TGase2 Polyclonal Antibody detects endogenous levels of TGase2 protein.
Storage and Stability:	-20°C/1 year