

Antibodies Product Datasheet

Product Name:	FUT8 Antibody
Catalog Number:	C48502
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FUT8
Concentration:	1 mg/ml
Gene Name:	FUT8
Protein Name:	Alpha-(1,6)-fucosyltransferase (alpha1-6FucT) (EC 2.4.1.68) (Fucosyltransferase 8) (GDP-L-Fuc:N-acetyl-beta-D-glucosaminide alpha1,6-fucosyltransferase) (GDP-fucose--glycoprotein fucosyltransferase) (Glycoprotein 6-alpha-L-fucosyltransferase)
Other Name:	
Human Swiss Prot No.:	Q9BYC5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FUT8 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year