

Antibodies Product Datasheet

Product Name:	HER4/ErbB4 (Phospho-Tyr984) Antibody
Catalog Number:	A9524
Reactivity:	Human, Rat, Mouse
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized phospho peptide around human HER4 and ErbB4 (Tyr984)
Concentration:	1 mg/ml
Gene Name:	ERBB4 HER4
Protein Name:	HER4/ErbB4 (Tyr984)
Other Name:	Receptor tyrosine-protein kinase erbB-4 (EC 2.7.10.1) (Proto-oncogene-like protein c-ErbB-4) (Tyrosine kinase-type cell surface receptor HER4) (p180erbB4) [Cleaved into: ERBB4 intracellular domain (4ICD) (E4ICD) (s80HER4)]
Human Swiss Prot No.:	Q15303
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human HER4/ErbB4 (phospho-Tyr984)
Storage and Stability:	-20°C/1 year