

Antibodies Product Datasheet

Product Name:	ELNE Antibody
Catalog Number:	C43606
Reactivity:	Human, Mouse, Bovine
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	ELANE ELA2
Protein Name:	Neutrophil elastase (EC 3.4.21.37) (Bone marrow serine protease) (Elastase-2) (Human leukocyte elastase) (HLE) (Medullasin) (PMN elastase)
Other Name:	
Human Swiss Prot No.:	P08246
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	ELNE Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year