

Antibodies Product Datasheet

Product Name:	HDAC6 (Phospho-Ser22) Antibody
Catalog Number:	A0941
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HDAC6 around the phosphorylation site of Ser22. AA range:7-56
Concentration:	1 mg/ml
Gene Name:	HDAC6
Protein Name:	Histone deacetylase 6
Other Name:	HDAC6; KIAA0901; JM21; Histone deacetylase 6; HD6
Human Swiss Prot No.:	Q9UBN7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-HDAC6 (S22) Polyclonal Antibody detects endogenous levels of HDAC6 protein only when phosphorylated at S22.
Storage and Stability:	-20°C/1 year