

## Antibodies Product Datasheet

<b>Product Name:</b>	Anti-Human GCDFP-15 (ABT075) Mouse Monoclonal Antibody
<b>Catalog Number:</b>	V0348
<b>Clone ID:</b>	ABT075
<b>Reactivity:</b>	Human
<b>Applications:</b>	IHC
<b>Application Dilutions</b>	IHC 1:100-200
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>UniProt:</b>	PIP
<b>Gene Name:</b>	P12273
<b>Other Name:</b>	BRST-2, BRST2, GCDFP-15, GCDGP15, gp17, GPIP4, Gross cystic disease fluid protein 15, PIP, PIP_HUMAN, Prolactin induced protein, Prolactin inducible protein precursor, Prolactin-induced protein, Prolactin-inducible protein, SABP, Secretory actin binding protein, Secretory actin-binding protein
<b>Isotype</b>	IgG1, Kappa
<b>Localization:</b>	Breast
<b>Positive Tissue Control:</b>	Cytoplasmic
<b>Immunogen</b>	Synthetic peptide
<b>Purification</b>	Immunogen affinity purified
<b>Formulation:</b>	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
<b>Specificity:</b>	The antibody can specifically recognize human GCDFP-15 protein.
<b>Storage and Stability:</b>	-20°C/1 year