

Antibodies Product Datasheet

Product Name:	Troponin T-C Antibody
Catalog Number:	C30058
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human TNNT2. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	TNNT2
Protein Name:	Troponin T cardiac muscle
Other Name:	TNNT2; Troponin T, cardiac muscle; TnTc; Cardiac muscle troponin T; cTnT
Human Swiss Prot No.:	P45379
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Troponin T-C Polyclonal Antibody detects endogenous levels of Troponin T-C protein.
Storage and Stability:	-20°C/1 year