

---

## Antibodies Product Datasheet

<b>Product Name:</b>	Elk-1 Antibody
<b>Catalog Number:</b>	B7070
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	IHC-p, IF(paraffin section), WB, ELISA
<b>Dilutions:</b>	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human Elk1. AA range:379-428
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	ELK1
<b>Protein Name:</b>	ETS domain-containing protein Elk-1
<b>Other Name:</b>	ELK1; ETS domain-containing protein Elk-1
<b>Human Swiss Prot No.:</b>	P19419
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Elk-1 Polyclonal Antibody detects endogenous levels of Elk-1 protein.
<b>Storage and Stability:</b>	-20°C/1 year