

## Antibodies Product Datasheet

<b>Product Name:</b>	SETD2 Antibody
<b>Catalog Number:</b>	C48811
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human SETD2
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	SETD2 HIF1 HYPB KIAA1732 KMT3A SET2 HSPC069
<b>Protein Name:</b>	Histone-lysine N-methyltransferase SETD2 (EC 2.1.1.43) (HIF-1) (Huntingtin yeast partner B) (Huntingtin-interacting protein 1) (HIP-1) (Huntingtin-interacting protein B) (Lysine N-methyltransferase 3A) (SET domain-containing protein 2) (hSET2) (p231HBP)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q9BYW2
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of SETD2 at Human
<b>Storage and Stability:</b>	Store at -20° C/1 year