

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Dlx-4 Antibody

Catalog Number: C11450

Reactivity: Human, Rat, Mouse

Applications: WB, IHC-p, IF/ICC, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in

other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human DLX4. AA range:61-110

Concentration: 1 mg/ml

Gene Name: DLX4

Protein Name: Homeobox protein DLX-4

Other Name: DLX4; BP1; DLX7; DLX8; DLX9; Homeobox protein DLX-4; Beta protein

1; Homeobox protein DLX-7; Homeobox protein DLX-8

Human Swiss Prot No.: Q92988

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Dlx-4 Polyclonal Antibody detects endogenous levels of Dlx-4 protein.

Storage and Stability: -20°C/1 year