
Antibodies Product Datasheet

Product Name:	BORA (Phospho-Ser497) Antibody
Catalog Number:	A8304
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human BORA (Phospho-Ser497)
Concentration:	1 mg/ml
Gene Name:	BORA C13orf34
Protein Name:	BORA (Phospho-Ser497)
Other Name:	Protein aurora borealis (HsBora)
Human Swiss Prot No.:	Q6PGQ7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.209% sodium azide.
Specificity:	This antibody detects endogenous levels of BORA (Phospho-Ser497) at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year