

Antibodies Product Datasheet

Product Name:	Cdk11A/B Antibody
Catalog Number:	C11909
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Cdk11A/B . at AA range: 190-270
Concentration:	1 mg/ml
Gene Name:	CDK11A/CDK11B
Protein Name:	Cell division protein kinase 11A/B
Other Name:	CDK11A; CDC2L2; CDC2L3; PITSLREB; Cyclin-dependent kinase 11A; Cell division cycle 2-like protein kinase 2; Cell division protein kinase 11A; Galactosyltransferase-associated protein kinase p58/GTA; PITSLRE serine/threonine-protein kinase C
Human Swiss Prot No.:	Q9UQ88/P21127
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdk11A/B Polyclonal Antibody detects endogenous levels of Cdk11A/B protein.
Storage and Stability:	-20°C/1 year