

## Antibodies Product Datasheet

<b>Product Name:</b>	His-tag Antibody
<b>Catalog Number:</b>	C40007
<b>Reactivity:</b>	Species independent
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	Western Blot: 1/1000 - 1/3000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	6X His synthetic peptide conjugated to KLH.
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	His-Tag
<b>Protein Name:</b>	His Tag
<b>Other Name:</b>	6 His epitope tag; Hexa His tag; HHHHHH epitope tag; HHHHHH tag; His tag
<b>Human Swiss Prot No.:</b>	
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.05% sodium azide.
<b>Specificity:</b>	His-tag Polyclonal Antibody detects His-tagged recombinant proteins or His-tagged proteins overexpressed in cells.
<b>Storage and Stability:</b>	-20°C/1 year