

Antibodies Product Datasheet

Product Name:	GRK 3 Antibody
Catalog Number:	C10560
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRK3. AA range:361-410
Concentration:	1 mg/ml
Gene Name:	ADRBK2
Protein Name:	Beta-adrenergic receptor kinase 2
Other Name:	ADRBK2; BARK2; GRK3; Beta-adrenergic receptor kinase 2; Beta-ARK-2; G-protein-coupled receptor kinase 3
Human Swiss Prot No.:	P35626
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GRK 3 Polyclonal Antibody detects endogenous levels of GRK 3 protein.
Storage and Stability:	-20°C/1 year