

Antibodies Product Datasheet

Product Name:	S4A7 Antibody
Catalog Number:	C46121
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from part region of human protein
Concentration:	1 mg/ml
Gene Name:	SLC4A7 BT NBC2 NBC2B NBC3 SBC2 SLC4A6
Protein Name:	Sodium bicarbonate cotransporter 3 (Sodium bicarbonate cotransporter 2) (Sodium bicarbonate cotransporter 2b) (Bicarbonate transporter) (Solute carrier family 4 member 7)
Other Name:	
Human Swiss Prot No.:	Q9Y6M7
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	S4A7 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year