

Antibodies Product Datasheet

Product Name:	IL-20R α Antibody
Catalog Number:	R12-2207
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IL-20R α . AA range:173-222
Concentration:	1 mg/ml
Gene Name:	IL20RA
Protein Name:	Interleukin-20 receptor subunit alpha
Other Name:	IL20RA; Interleukin-20 receptor subunit alpha; IL-20 receptor subunit alpha; IL-20R-alpha; IL-20RA; Cytokine receptor class-II member 8; Cytokine receptor family 2 member 8; CRF2-8; IL-20R1; ZcytoR7
Human Swiss Prot No.:	Q9UHF4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-20R α Polyclonal Antibody detects endogenous levels of IL-20R α protein.
Storage and Stability:	-20°C/1 year