

Antibodies Product Datasheet

Product Name:	eIF1AY Antibody
Catalog Number:	C15692
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EIF1AY. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	EIF1AY
Protein Name:	Eukaryotic translation initiation factor 1A Y-chromosomal
Other Name:	EIF1AY; Eukaryotic translation initiation factor 1A; Y-chromosomal; eIF-1A Y isoform; Eukaryotic translation initiation factor 4C; eIF-4C
Human Swiss Prot No.:	O14602
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	eIF1AY Polyclonal Antibody detects endogenous levels of eIF1AY protein.
Storage and Stability:	-20°C/1 year