

Antibodies Product Datasheet

Product Name:	Cleaved-Integrin α 5 HC (F42) Antibody
Catalog Number:	L0277
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ITGA5. AA range:23-72
Concentration:	1 mg/ml
Gene Name:	ITGA5
Protein Name:	Integrin alpha-5
Other Name:	ITGA5; FNRA; Integrin alpha-5; CD49 antigen-like family member E; Fibronectin receptor subunit alpha; Integrin alpha-F; VLA-5; CD antigen CD49e
Human Swiss Prot No.:	P08648
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Integrin α 5 HC (F42) Polyclonal Antibody detects endogenous levels of fragment of activated Integrin α 5 HC protein resulting from cleavage adjacent to F42.
Storage and Stability:	-20°C/1 year