

Antibodies Product Datasheet

Product Name:	SENP2 mouse mAb(PRT2320)
Catalog Number:	N1750
Reactivity:	Human, Mouse,Rat
Applications:	WB
Dilutions:	WB 1:1000-2000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human SENP2
Concentration:	1 mg/ml
Gene Name:	SENP2 KIAA1331
Protein Name:	SENP2
Other Name:	Sentrin-specific protease 2 (EC 3.4.22.68) (Axa2) (SMT3-specific isopeptidase 2) (Smt3ip2) (Sentrin/SUMO-specific protease SENP2)
Human Swiss Prot No.:	Q9HC62
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human SENP2
Storage and Stability:	-20°C/1 year