

Antibodies Product Datasheet

Product Name:	Anti-Cytokeratin 20 [ABT-Ck20]
Catalog Number:	V0052
Clone ID:	ABT-Ck20
Reactivity:	Human
Applications:	IHC
Application Dilutions	IHC 1:100-200
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	P35900
Gene Name:	KRT20
Other Name:	Keratin, type I cytoskeletal 20; Cytokeratin-20; CK-20; K20
Isotype	Mouse IgG
Localization:	N/A
Positive Tissue Control:	N/A
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	Anti-Cytokeratin 20 [ABT-Ck20] detects endogenous levels of total Cytokeratin 20 protein.
Storage and Stability:	-20°C/1 year