

## Antibodies Product Datasheet

<b>Product Name:</b>	β-Catenin Antibody
<b>Catalog Number:</b>	C40034
<b>Reactivity:</b>	Mouse, Rat
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB: 1:2000-4000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of Catenin beta-1
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	CTNNB1
<b>Protein Name:</b>	Catenin beta-1
<b>Other Name:</b>	CTNNB1; CTNNB; OK/SW-cl.35; Catenin beta-1; Beta-catenin; YM3345
<b>Human Swiss Prot No.:</b>	P35222
<b>Formulation:</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous β-Catenin protein.
<b>Storage and Stability:</b>	-20°C/1 year