

## Antibodies Product Datasheet

<b>Product Name:</b>	MGAT4C Antibody
<b>Catalog Number:</b>	C49055
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB
<b>Dilutions:</b>	WB 1:500-2000
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal, Rabbit,IgG
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human MGAT4C
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	MGAT4C
<b>Protein Name:</b>	Alpha-1,3-mannosyl-glycoprotein 4-beta-N-acetylglucosaminyltransferase C (EC 2.4.1.145) (N-acetylglucosaminyltransferase IV homolog) (hGnT-IV-H) (N-glycosyl-oligosaccharide-glycoprotein N-acetylglucosaminyltransferase IVc) (GlcNAc-T IVc) (GnT-IVc) (N-acet
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	Q9UBM8
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of MGAT4C at Human, Mouse
<b>Storage and Stability:</b>	Store at -20° C/1 year