
Antibodies Product Datasheet

Product Name:	Rap 2C Antibody
Catalog Number:	C18282
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAP2C. AA range:97-146
Concentration:	1 mg/ml
Gene Name:	RAP2C
Protein Name:	Ras-related protein Rap-2c
Other Name:	RAP2C; Ras-related protein Rap-2c
Human Swiss Prot No.:	Q9Y3L5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rap 2C Polyclonal Antibody detects endogenous levels of Rap 2C protein.
Storage and Stability:	-20°C/1 year