

Antibodies Product Datasheet

Product Name:	TLR2 Antibody
Catalog Number:	C47553
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section)
Dilutions:	IHC: 1:200-500
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of TLR2
Concentration:	
Gene Name:	TLR2
Protein Name:	Toll-like receptor 2
Other Name:	TLR2; TIL4; Toll-like receptor 2; Toll/interleukin-1 receptor-like protein 4; CD282; YM3391
Human Swiss Prot No.:	O60603
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous TLR2 protein.
Storage and Stability:	-20°C/1 year