
Antibodies Product Datasheet

Product Name:	ATPAF2 Antibody
Catalog Number:	C14600
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATPAF2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	ATPAF2
Protein Name:	ATP synthase mitochondrial F1 complex assembly factor 2
Other Name:	ATPAF2; ATP12; LP3663; ATP synthase mitochondrial F1 complex assembly factor 2; ATP12 homolog
Human Swiss Prot No.:	Q8N5M1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATPAF2 Polyclonal Antibody detects endogenous levels of ATPAF2 protein.
Storage and Stability:	-20°C/1 year