

## **Assay Biotechnology Company**

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

## **Antibodies Product Datasheet**

**Product Name:** ABHD7 Antibody

Catalog Number: C14219

Reactivity: Human, Mouse

**Applications:** WB, ELISA

Dilutions: Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other

applications.

**Modification:** Unmodified/Total

**Source:** Polyclonal Rabbit

**Purification:** The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

**Immunogen:** The antiserum was produced against synthesized peptide derived from

human ABHD7. AA range:101-150

Concentration: 1 mg/ml

Gene Name: EPHX4

**Protein Name:** Epoxide hydrolase 4

Other Name: EPHX4; ABHD7; EPHXRP; Epoxide hydrolase 4; Abhydrolase domain-

containing protein 7; Epoxide hydrolase-related protein

Human Swiss Prot No.: Q8IUS5

**Formulation:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

**Specificity:** ABHD7 Polyclonal Antibody detects endogenous levels of ABHD7

protein.

Storage and Stability: -20°C/1 year