

Antibodies Product Datasheet

Product Name:	CREB Antibody
Catalog Number:	C47497
Reactivity:	Human
Applications:	WB
Dilutions:	WB: 1:1000-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	Recombinant Protein of Cyclic AMP-responsive element-binding protein 1
Concentration:	
Gene Name:	CREB1
Protein Name:	Cyclic AMP-responsive element-binding protein 1
Other Name:	CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1; YM3239
Human Swiss Prot No.:	P16220
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous CREB protein.
Storage and Stability:	-20°C/1 year