

## Antibodies Product Datasheet

<b>Product Name:</b>	3pK Antibody
<b>Catalog Number:</b>	C11133
<b>Reactivity:</b>	Human, Mouse, Rat
<b>Applications:</b>	WB, IHC-p, IF/ICC, ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/Total
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human MAPK3. AA range:301-350
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	MAPKAPK3
<b>Protein Name:</b>	MAP kinase-activated protein kinase 3
<b>Other Name:</b>	MAPKAPK3; MAP kinase-activated protein kinase 3; MAPK-activated protein kinase 3; MAPKAP kinase 3; MAPKAP-K3; MAPKAPK-3; MK-3; Chromosome 3p kinase; 3pK
<b>Human Swiss Prot No.:</b>	Q16644
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	3pK Polyclonal Antibody detects endogenous levels of 3pK protein.
<b>Storage and Stability:</b>	-20°C/1 year