
Antibodies Product Datasheet

Product Name:	Aldolase A Antibody
Catalog Number:	C10268
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ALDOA. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	ALDOA
Protein Name:	Fructose-bisphosphate aldolase A
Other Name:	ALDOA; ALDA; Fructose-bisphosphate aldolase A; Lung cancer antigen NY-LU-1; Muscle-type aldolase
Human Swiss Prot No.:	P04075
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Aldolase A Polyclonal Antibody detects endogenous levels of Aldolase A protein.
Storage and Stability:	-20°C/1 year