

Antibodies Product Datasheet

Product Name:	ERCC4 Antibody
Catalog Number:	C0394
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human XPF. AA range:801-850
Concentration:	1 mg/ml
Gene Name:	ERCC4
Protein Name:	DNA repair endonuclease XPF
Other Name:	ERCC4; ERCC11; XPF; DNA repair endonuclease XPF; DNA excision repair protein ERCC-4; DNA repair protein complementing XP-F cells; Xeroderma pigmentosum group F-complementing protein
Human Swiss Prot No.:	Q92889
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ERCC4 Polyclonal Antibody detects endogenous levels of ERCC4 protein.
Storage and Stability:	-20°C/1 year