

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: CD95 Antibody

Catalog Number: C0355

Reactivity: Human, Rat, Mouse

Applications: WB, IHC-p, IF(paraffin section), ELISA

Dilutions: Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/20000. Not yet tested in other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human FAS. AA range:281-330

Concentration: 1 mg/ml

Gene Name: FAS

Protein Name: Tumor necrosis factor receptor superfamily member 6

Other Name: FAS; APT1; FAS1; TNFRSF6; Tumor necrosis factor receptor

superfamily member 6; Apo-1 antigen; Apoptosis-mediating surface

antigen FAS; FASLG receptor; CD antigen CD95

Human Swiss Prot No.: P25445

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: CD95 Polyclonal Antibody detects endogenous levels of CD95 protein.

Storage and Stability: -20°C/1 year