

Antibodies Product Datasheet

Product Name:	Rit1 Antibody
Catalog Number:	C0317
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RIT1. AA range:170-219
Concentration:	1 mg/ml
Gene Name:	RIT1
Protein Name:	GTP-binding protein Rit1
Other Name:	RIT1; RIBB; RIT; ROC1; GTP-binding protein Rit1; Ras-like protein expressed in many tissues; Ras-like without CAAX protein 1
Human Swiss Prot No.:	Q92963
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rit1 Polyclonal Antibody detects endogenous levels of Rit1 protein.
Storage and Stability:	-20°C/1 year