

Antibodies Product Datasheet

Product Name:	Oct-2 Antibody
Catalog Number:	C0282
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human OCT2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	POU2F2
Protein Name:	POU domain class 2 transcription factor 2
Other Name:	POU2F2; OCT2; OTF2; POU domain; class 2, transcription factor 2; Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2; Octamer-binding protein 2; Oct-2; Octamer-binding transcription factor 2; OTF-2
Human Swiss Prot No.:	P09086
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Oct-2 Polyclonal Antibody detects endogenous levels of 40818 protein.
Storage and Stability:	-20°C/1 year