

Assay Biotechnology Company

47787 Fremont Blvd. Fremont, CA 94538 Phone: (408) 747 0185 Fax: (408 747 0145 www.assaybiotechnlogy.com

Antibodies Product Datasheet

Product Name: Ki-67 Antibody

Catalog Number: C0250

Reactivity: Human, Rat, Mouse

Applications: WB, IHC-p, IF/ICC, ELISA

Dilutions: Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 -

1/1000. ELISA: 1/20000. Not yet tested in other applications.

Modification: Unmodified/Total

Source: Polyclonal Rabbit

Purification: The antibody was affinity-purified from rabbit antiserum by affinity-

chromatography using epitope-specific immunogen.

Immunogen: The antiserum was produced against synthesized peptide derived from

human Ki67. AA range:3207-3256

Concentration: 1 mg/ml

Gene Name: MKI67

Protein Name: Antigen KI-67

Other Name: MKI67; Antigen KI-67

Human Swiss Prot No.: P46013

Formulation: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium

azide.

Specificity: Ki-67 Polyclonal Antibody detects endogenous levels of Ki-67 protein.

Storage and Stability: -20°C/1 year