

Antibodies Product Datasheet

Product Name:	HSP90A Antibody
Catalog Number:	C0234
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HSP90A. AA range:683-732
Concentration:	1 mg/ml
Gene Name:	HSP90AA1
Protein Name:	Heat shock protein HSP 90-alpha
Other Name:	HSP90AA1; HSP90A; HSPC1; HSPCA; Heat shock protein HSP 90-alpha; Heat shock 86 kDa; HSP 86; HSP86; Renal carcinoma antigen NY-REN-38
Human Swiss Prot No.:	P07900
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HSP90A Polyclonal Antibody detects endogenous levels of HSP90A protein.
Storage and Stability:	-20°C/1 year