

Antibodies Product Datasheet

Product Name:	YB-1 Antibody
Catalog Number:	B8481
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300.Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity- chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human YBX1. AA range:68-117
Concentration:	1 mg/ml
Gene Name:	YBX1
Protein Name:	Nuclease-sensitive element-binding protein 1
Other Name: Human Swiss Prot No.:	YBX1; NSEP1; YB1; Nuclease-sensitive element-binding protein 1; CCAAT-binding transcription factor I subunit A; CBF-A; DNA-binding protein B; DBPB; Enhancer factor I subunit A; EFI-A; Y-box transcription factor; Y-box-binding protein 1; YB- P67809
Human Swiss Flot No	F 07 609
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	YB-1 Polyclonal Antibody detects endogenous levels of YB-1 protein.
Storage and Stability:	-20°C/1 year