

Antibodies Product Datasheet

Product Name:	ER β Antibody
Catalog Number:	B8280
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ER β . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	ESR2
Protein Name:	Estrogen receptor beta
Other Name:	ESR2; ESTRB; NR3A2; Estrogen receptor beta; ER-beta; Nuclear receptor subfamily 3 group A member 2
Human Swiss Prot No.:	Q92731
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ER β Polyclonal Antibody detects endogenous levels of ER β protein.
Storage and Stability:	-20°C/1 year